

Signal-Seeker PTM Detection Kits

15% OFF



Select Signal-Seeker Products

Comprehensive kits

Minimal optimization time

Maximizes PTM detection

New Acetyl Lysine Affinity-Beads and Kits

VALID UNTIL



2025

EUROMEDEX

24 rue des Tuileries, 67460 Souffelweyersheim

Téléphone : +33 (0)3 88 18 07 22

 Visiter la boutique  Suivez-nous sur LinkedIn



Popular Signal-Seeker Tools

Product Summary Highlights


Acetyl-Lysine antibody (CS-AAC04-Beads)

Simplified acetylation detection with Cytoskeleton's Acetyl Lysine Affinity Beads-AAC04

Fig Legend: AAC03 & AAC02 acetyl lysine affinity beads (50µl bead slurry) were used to IP acetylated proteins from Cos-7 cells either treated (+) or untreated (-) with deacetylase inhibitors.

BUY

Learn more at: <https://www.cytoskeleton.com/aac03-beads>



	AAC03	AAC02
Acetyl-EGFR Signal ratio (+/-)	4.0x	1.4x
Acetyl RhoGDI Signal ratio (+/-)	2.8x	1.0x


SUMO1 affinity beads (CS-ASM11-beads)

The only SUMO1 affinity beads commercially available with precise specificity for SUMO1 modified proteins

BUY

Investigate endogenous SUMO1 modifications for any target protein

Learn more at: <https://www.cytoskeleton.com/sumo1-asm11-beads>



Signal-Seeker Ubiquitin Detection Kit (CS-BK161)

The original, comprehensive ubiquitin detection kit: contains all of the essential components needed for effective ubiquitination detection

Perfect for scientists that are new to the Ubiquitin field

BUY

Learn more at: <https://www.cytoskeleton.com/bk161>

Essential Kit Components:

- IP Affinity Beads: (Cat. # UBA01-beads)
- IP control Beads: (Cat. # CUB02-beads)
- Universal Lysis System
- Key Inhibitors
- Wash Reagents
- Elution Reagents
- Spin Columns
- Protein Quantitation System
- Chemiluminescence

EUROMEDEX

24 rue des Tuileries, 67460 Souffelweyersheim
Téléphone : +33 (0)3 88 18 07 22

Visiter la boutique Suivez-nous sur LinkedIn

Ref Number	Product	Size	List Price € HT	Offer Price € HT
CS-MXA95	Actin Protein (MICAL-Oxidized): Rabbit Skeletal Muscle (Rabbit Muscle, > 95% Pure)	2 x 250 µg	469,00	398,65
CS-MXA95-XL	Actin Protein (MICAL-Oxidized): Rabbit Skeletal Muscle (Rabbit Muscle, > 95% Pure)	1 x 1 mg	756,00	642,60
CS-AAC03-S	Anti-Acetyl Lysine Mouse Monoclonal Antibody (Trial size MAb clone 19C4B2.1)	1 x 25 µl	150,00	127,50
CS-AAC02-S	Anti-Acetyl Lysine Mouse Monoclonal Antibody (Trial size MAb clone 7B5A1)	1 x 25 µl	150,00	127,50
CS-AAC02-beads	Anti-Acetyl Lysine Mouse Monoclonal Antibody Agarose Conjugate. MAb clone 7B5A1	4 x 500 µl	848,00	720,80
CS-AAC03-beads	Anti-Acetyl Lysine Mouse Monoclonal Antibody Agarose Conjugates MAb clone 19C4B2.1	4 x 500 µl	848,00	720,80
CS-AAC04-beads	Anti-Acetyl Lysine Mouse Monoclonal Antibody Agarose Conjugates MAb clones 19C4B2.1 and 7B5A1	4 x 500 µl	848,00	720,80
CS-AAC03-HRP	Anti-Acetyl Lysine Mouse Monoclonal Antibody HRP Conjugates (MAb Clone 19C4B2.1)	1 x 100 µl	606,00	515,10
CS-AAC03-HRP-S	Anti-Acetyl Lysine Mouse Monoclonal Antibody HRP Conjugates (Trial size MAb clone 19C4B2.1)	1 x 25 µl	209,00	177,65
CS-AAC03	Anti-Acetyl Lysine Mouse Monoclonal Antibody MAb clone 19C4B2.1	2 x 100 µl	541,00	459,85
CS-AAC01	Anti-Acetyl Lysine Mouse Monoclonal Antibody MAb clone 3C6.08.20	2 x 100 µl	541,00	459,85
CS-AAC01-S	Anti-Acetyl Lysine Mouse Monoclonal Antibody MAB clone 3C6.08.20	1 x 25 µl	150,00	127,50
CS-AAC02	Anti-Acetyl Lysine Mouse Monoclonal Antibody MAb clone 7B5A1	2 x 100 µl	541,00	459,85
CS-APY03-Beads	Anti-Phosphotyrosine Affinity Beads	4 x 330 µl	848,00	720,80
CS-APY03	Anti-Phosphotyrosine Mouse Monoclonal Antibody	2 x 100 µl	541,00	459,85
CS-APY03-S	Anti-Phosphotyrosine Mouse Monoclonal Antibody (Trial Size)	1 x 25 µl	150,00	127,50
CS-APY03-HRP	Anti-Phosphotyrosine-HRP Mouse Monoclonal Antibody	1 x 100 µl	606,00	515,10
CS-APY03-HRP-S	Anti-Phosphotyrosine-HRP Mouse Monoclonal Antibody (Trial Size)	1 x 25 µl	202,00	171,70
CS-ASM01	Anti-SUMO-1 Mouse Monoclonal Antibody	1 x 100 µl	541,00	459,85
CS-ASM01-S	Anti-SUMO-1 Mouse Monoclonal Antibody (Trial Size)	1 x 25 µl	150,00	127,50
CS-ASM23	Anti-SUMO-2/3 Mouse Monoclonal Antibody	2 x 100 µl	541,00	459,85
CS-ASM23-HRP	Anti-SUMO-2/3 Mouse Monoclonal Antibody	1 x 100 µl	613,00	521,05
CS-ASM24	Anti-SUMO-2/3 Mouse Monoclonal Antibody	2 x 200 µl	541,00	459,85
CS-ASM23-HRP-S	Anti-SUMO-2/3 Mouse Monoclonal Antibody (Trial Size)	1 x 25 µl	196,00	166,60
CS-ASM23-S	Anti-SUMO-2/3 Mouse Monoclonal Antibody (Trial Size)	1 x 25 µl	150,00	127,50
CS-ASM24-S	Anti-SUMO-2/3 Mouse Monoclonal Antibody (Trial Size)	1 x 150 µl	196,00	166,60
CS-AUB01	Anti-Ubiquitin Mouse Monoclonal Antibody	2 x 100 µl	541,00	459,85
CS-AUB01-XL	Anti-Ubiquitin Mouse Monoclonal Antibody	4 x 500 µl	1154,00	980,90
CS-AUB01-S	Anti-Ubiquitin Mouse Monoclonal Antibody (Trial size)	1 x 25 µl	150,00	127,50
CS-AUB01-HRP	Anti-Ubiquitin-HRP Antibody	1 x 100 µl	613,00	521,05
CS-AUB01-HRP-S	Anti-Ubiquitin-HRP Antibody	1 x 25 µl	196,00	166,60
CS-BLR03	BlastR Lysis/Dilution Kit	50 lysates	202,00	171,70
CS-BLR02	BlastR Rapid Filtration Kit	50 purifications	202,00	171,70
CS-BLR01	BlastR Rapid Lysate Prep Kit	50 purifications	306,00	260,10
CS-AACX1	Chemically Acetylated Bovine Serum Albumin (BSA)	2 x 500 µg	78,00	66,30
CS-CIG02-beads	Control for Acetylation Immunoprecipitation IgG beads	10 assays	150,00	127,50
CS-CIG01-beads	Control for Immunoprecipitation IgG beads	10 assays	150,00	127,50
CS-CIG03-beads	Control for SUMO1 Immunoprecipitation IgG beads	10 assays	150,00	127,50
CS-CUB02B-beads	Control for Ubiquitin Affinity Beads	10 assays	146,00	124,10

CS-MIC01	MICAL-1 Protein Redox-CH domains (6xHis-tag)	2 x 50 µg	495,00	420,75
CS-MIC01-XL	MICAL-1 Protein Redox-CH domains (6xHis-tag)	1 x 1 mg	1878,00	1596,30
CS-MB201	MsrB2 Protein aa: 42-182 (6xHis-tag)	2 x 50 µg	443,00	376,55
CS-MB201-XL	MsrB2 Protein aa: 42-182 (6xHis-tag)	1 x 1 mg	1878,00	1596,30
CS-MPAX1	Pyrene Labeled Actin Protein (MICAL-Oxidized): Rabbit Skeletal Muscle (Rabbit Muscle, > 95% Pure)	2 x 250 µg	548,00	465,80
CS-MPAX1-XL	Pyrene Labeled Actin Protein (MICAL-Oxidized): Rabbit Skeletal Muscle (Rabbit Muscle, > 95% Pure)	1 x 1 mg	887,00	753,95
CS-BK164	Signal-Seeker® Acetyl-Lysine Detection Kit - AAC03-Beads (30 assays)	KIT 30 assays	1030,00	875,50
CS-BK163-S	Signal-Seeker® Acetyl-Lysine Detection Kit (10 assays)	KIT 10 assays	520,00	442,00
CS-BK163	Signal-Seeker® Acetyl-Lysine Detection Kit (30 assays)	KIT 30 assays	1030,00	875,50
CS-BK160-S	Signal-Seeker® Phosphotyrosine Detection Kit (10 assays)	KIT 10 assays	520,00	442,00
CS-BK160	Signal-Seeker® Phosphotyrosine Detection Kit (30 assays)	KIT 30 assays	1030,00	875,50
CS-BK162-S	Signal-Seeker® SUMOylation Detection Kit (10 assays)	KIT 10 assays	520,00	442,00
CS-BK162	Signal-Seeker® SUMOylation Detection Kit (30 assays)	KIT 30 assays	1030,00	875,50
CS-BK165-S	Signal-Seeker® SUMOylation-1 Detection Kit (10 assays)	KIT 10 assays	520,00	442,00
CS-BK165	Signal-Seeker® SUMOylation-1 Detection Kit (30 assays)	KIT 30 assays	1030,00	875,50
CS-BK161-S	Signal-Seeker® Ubiquitin Detection Kit (10 assays)	KIT 10 assays	520,00	442,00
CS-BK161	Signal-Seeker® Ubiquitin Detection Kit (30 assays)	KIT 30 assays	1030,00	875,50
CS-ASM24-Beads	SUMO-2/3 Affinity Beads	2 x 400 µl	848,00	720,80
CS-ASM11-Beads	SUMOylation-1 Affinity Beads	3 x 400 µl	848,00	720,80
CS-AACX2	Tricostatin A (TSA) Treated And Untreated Cell Lysates	100 µg treated + 100 µg untreated	150,00	127,50
CS-UBA01B-beads	Ubiquitin Affinity Beads	40 rxs	750,00	637,50